

Professional opinion about the HEINE® ML 4 LED HeadLight



Dr. Lisa Berchthold

The use of headlights has become standard for the dental care of horses and is meanwhile indispensable for VETs. Bad light conditions for treating veterinarians in cowsheds or barns are however still the rule rather than the exception. Using a headlight greatly improves working conditions in the equine and large animal practice – and can also be used for surgical wound treatments or outpatient surgery in sheds. A cesarean section in a cow is often performed at night where a good headlight improves the work of the veterinarian and speeds up surgery. The shadow-free illumination of the head worn light also provides optimal support when performing surgical claw treatments or taking care of teat injuries.

The HEINE ML 4 LED HeadLight is particularly well suited for veterinarians in the equine and large animal practice and is a useful accessory for every practice car. The user benefits in particular from the cable-free rechargeable battery, which is integrated into the headband. It is therefore no obstacle to winter cloths, surgical or delivery gowns. The rechargeable battery has a continuous operation time of 3.5 hours.

The height and width of the headband is individually adjustable at two points and provides perfect and



secure comfort of wear. This combined with the integrated rechargeable battery allows the vet to move freely in the shed and in between animals.

The brightness of the ML 4 LED HeadLight is continuously adjustable and the precision multiple lens optical system consisting of four high-quality lenses provides an absolute homogeneous illumination of the relevant examination, treatment and operation areas.

Dr. Lisa Berchthold is a certified veterinarian since 2006 and works in a mixed veterinary practice in the district of Weilheim, Germany since 2009.



HEINE ML 4 LED HeadLight

